


RACKS, BOXES & CRYOVIALS



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Россия +7(495)268-04-70

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новыйбурск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47
Казахстан +7(7172)727-132

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35
Киргизия +996(312)96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

CRYO BOXES


Arctiko Cryo boxes are available in 2 versions. They are made of high quality and perfect for short- or long-term storage. Our standard 9x9 boxes fits perfectly in various rack systems and has a storage capacity of 81 samples. The new and innovative 10x10 box are perfectly matched together with the new and smart Arctiko Cryovials that together makes it possible to store 100 samples in only 1 box.

The 10x10 dividers provides 23% additional storage capacity compared to our conventional solution. The boxes are made of high quality durable

polypropylene and can be used in temperatures as low as -150°C^* . The boxes are autoclavable at $+121^{\circ}\text{C}$ for 20 minutes.

With features like graduated alpha-numeric coding, smart closing indication and high-quality durable resin, these boxes are made with the purpose of fast and easy handling.

* Internally tested to be able to withstand down to -150°C .



PRODUCT OVERVIEW

The clear lid is connected to the box with special lock 'click on system' to prevent accidental opening. The upper corners of the lid and base are cut with an angle for faster identification and handling.

Optically clear lid and colored base with openings for ventilation.

Numerical coding for each tube enables quick and easy sample identification.

Arctiko 10x10 Cryo boxes fits perfectly with Arctiko cryo vials: 0.5, 1.0, 2.0 and 4.0.

A1	A2	A3	A4	A5	A6	A7	A8	A9	A10
B1	B2	B3	B4	B5	B6	B7	B8	B9	B10
C1	C2	C3	C4	C5	C6	C7	C8	C9	C10
D1	D2	D3	D4	D5	D6	D7	D8	D9	D10
E1	E2	E3	E4	E5	E6	E7	E8	E9	E10
F1	F2	F3	F4	F5	F6	F7	F8	F9	F10
G1	G2	G3	G4	G5	G6	G7	G8	G9	G10
H1	H2	H3	H4	H5	H6	H7	H8	H9	H10
I1	I2	I3	I4	I5	I6	I7	I8	I9	I10
J1	J2	J3	J4	J5	J6	J7	J8	J9	J10

Example of a 81 vial (9x9 / left), and a 100 vial (10x10 / right) cryo box layout

PRODUCT FEATURES

- Connects to the base with a click to prevent accidental opening
- Upper edges cut at an angle to ensure proper and fast lid positioning (only 10x10 boxes)
- Alpha-numerical system on top and bottom
- Molded in for easy sample identification
- Autoclavable
- 4 openings for ventilation
- Star feet lock
- Enables one handed opening of vials without taking them out
- 130x130 DIMENSIONS FITS MOST RACK SYSTEMS



CRYOVIALS



23% MORE SAMPLES WITH THE USE OF ARCTIKO CRYOVIALS AND BOXES, DUE TO THE INNOVATIVE DESIGN.

Arctiko Cryovials is uniquely designed with external threads for safe short- and long-term storage. The cryovials are constructed of high quality special formulated polypropylene which makes the Arctiko cryovials highly resistant and able to withstand temperatures of -196 °C to 121 °C.

Arctiko's Cryovials combined with the innovative 10x10 Cryoboxes, makes it possible to store 100 pcs in only one box without any extra costs involved. Compared with the conventional box, you will get room for 19 more samples without having to compromise with the space.

PRODUCT FEATURES

- ▶ 100 pcs of Cryovials fits in the innovative 10x10 boxes
- ▶ More than 23% more samples
- ▶ External thread design
- ▶ Smart star-foot grip
- ▶ Uniquely designed lids



QUALITY CERTIFICATES

Arctiko's cryovials are designed and manufactured in Denmark ensuring the best possible quality.

- ▶ Quality tested
- ▶ Produced according to GMP guidelines in a cleanroom
- ▶ RNase-free, DNase-free & Pyrogen-free
- ▶ DNA, RNA & ATP-free
- ▶ Sterilized after ISO 11137-2:2013
- ▶ Sterility Assurance Level (SAL) by standard products 10^{-6} according to DIN EN ISO 11137
- ▶ 5 years from date of sterilization - As long as the packaging remains unopened and free
- ▶ Free from detectable heavy metals, BSE / TSE and cytotoxic substances from visible defect



SPECIFICATIONS

CRYOVIALS SPECIFICATIONS	CV05	CV10	CV20	CV40	CV50
Volume (ml)	0,5	1,0	2,0	4,0	5,0
Length (mm)	28,6	36,1	47,6	73,3	88,5
Diameter (mm)	12,08	12,08	12,08	12,08	12,08
Material	PP (Polypropylen)	PP (Polypropylen)	PP (Polypropylen)	PP (Polypropylen)	PP (Polypropylen)
Temperature range	-196°C to 121°C	-196°C to 121°C	-196°C to 121°C	-196°C to 121°C	-196°C to 121°C
Pressure tolerance	95 kPa (0.95 bar)	95 kPa (0.95 bar)	95 kPa (0.95 bar)	95 kPa (0.95 bar)	95 kPa (0.95 bar)
Max filling (Vapor LN ₂) -196 °C	0,4	0,8	1,3	3,3	4,3
Max filling -20 °C	0,5	0,9	1,7	3,5	4,5
Max filling -80 °C	0,4	0,9	1,5	3,5	4,5
Order quantity	1 BOX (10X100 PCS PR. BACK) = 1000 PCS	1 BOX (10X100 PCS PR. BACK) = 1000 PCS	1 BOX (10X100 PCS PR. BACK) = 1000 PCS	1 BOX (5X100 PCS PR. BACK) = 500 PCS	1 BOX (5X100 PCS PR. BACK) = 500 PCS


CRYO BOXES SPECIFICATIONS	PBN05	PBN10	PBN20	PBN40	PBN50	PB50	PB75	PB96
Item no.	DAI 3160	DAI 3161	DAI 3162	DAI 3163	DAI 3164	DAI 3100	DAI 3106	DAI 3112
Box size (Height in mm)	32	40	50	80	90	50	75	96
Dimensions (WxDxH in mm)	130x130x32	130x130x40	130x130x50	130x130x80	130x130x90	133x133x50	133x133x75	133x133x96
Total vial capacity	100 pcs.	100 pcs.	100 pcs.	100 pcs.	100 pcs.	81 pcs.	81 pcs.	81 pcs.
Order quantity	1 pcs (10x10 BOX)	1 pcs (10x10 BOX)	1 pcs (10x10 BOX)	1 pcs (10x10 BOX)	1 pcs (10x10 BOX)	1 pcs (9x9 BOX) 1 PCS (9x9 BOX)	1 pcs (9x9 BOX)	1 pcs (9x9 BOX)

RACKS

MAKE YOUR SCIENTIFIC LIFE EASIER.

Store your freezer storage boxes within organized freezer racks. Arctiko offers stainless steel racks to optimize your freezer space and makes it easy to find the items you desire. Don't waste your time sifting through unsorted samples.

We offer a large variety of racks for our freezers. These racks are made of high-quality stainless steel and help you organize the samples inside your freezer. Our range consists of racks for both chest freezers and upright freezers, either with drawers or with compartments in order to introduce the box from the side.



UPRIGHT DRAWER RACK

Strong stainless steel racks with drawers for easy and practical access to your boxes. The drawers slide smoothly and easily, which facilitates easy removal and insertion.

The upright drawer rack is equipped with a smart grip, which ensures easy handling and finding of necessary samples.

When using racks with drawers, only the relevant drawer with the freezer boxes of interest is pulled-out. All other boxes stay safely in the cold interior of the ULT freezer. Arctiko's UDR are designed for smooth and ergonomic hand-lining, as well as to give a great overview of samples. The upright drawer racks fit all boxes up to 133x133x37 mm.




1 Smart design that optimizes the freezer space

3 Open sides for optimal freezing of Samples

2 Smart grip for smooth handling

4 Smooth slide of the drawers



CHEST FRONT RACK

The chest front rack is designed for optimal suiting of chest freezers. Due to the design and choice of material, these unique racks are extremely solid and long-lasting stainless steel. These racks are easy to work with and give a good grip of the racks, thanks to the smart folding handle on top of the rack.

The box is removed from the front, thus ensuring a good overview. The handle in the top meanwhile ensures easy handling of the racks. The chest front rack fit all boxes up to 133x133x37 mm.



1
Folded handle in the top
for easy handling



2
Made of solid and
long-lasting stainless steel

3
Optimized sample space due
to the smart thin and low
weighted stainless steel

4
Fits perfectly for
chest freezers
(can handle up to 12 boxes)

UPRIGHT SIDE RACK

The upright side rack model features individual compartments for each box, meaning that the box is introduced from the side with the folding vertical handle.

The rack is smart in the event you need to access more samples at the time. Open on both sides of the rack for great airflow around the samples. The upright side racks fit all boxes up to 133x133x37 mm.



Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Россия +7(495)268-04-70

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Казахстан +7(7172)727-132

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Киргизия +996(312)96-26-47

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

ARCTIKO

true cooling specialists

